

Fluoroquinolone Immunoassay Gold-Labeled Rapid Test Strips



Cat. No. ART-013

Lot. No. (See product label)

Product Name

Fluoroquinolone Immunoassay Gold-Labeled Rapid Test Strips

Product Introduction

This product is intended for the qualitative and semi-quantitative detection of fluoroquinolone residues in tissue samples. The entire testing process takes approximately 5 to 8 minutes. The detection sensitivity reaches 50 µg/kg (ppb).

Applicable Samples

Tissue, Eggs, Honey.

Precautions

1. Please follow the operating procedures. Do not touch the test strip's display area during use.
2. If you discover expired, damaged, contaminated, or defective products upon purchase, please return them to the place of purchase for replacement.
3. This product is for single-use only. Use the test strip immediately after opening; do not reuse it.
4. Do not add samples to the sample well immediately after removing them from the refrigerator or after processing. Allow the sample to return to room temperature ($20\pm5^{\circ}\text{C}$) before testing.
5. Do not use water or other liquids as negative control samples.

Shelf Life

Shelf life is 12 months. The production date, expiration date, and lot number are indicated on the outer packaging.

Storage

Store the original packaging at 2-30°C in a cool, dry place away from light. Do not freeze or expose to direct sunlight.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY