

Porcine Toxoplasma Antibody Rapid Test Strip



Cat. No. ADT-017

Lot. No. (See product label)

Product Name

Porcine Toxoplasma Antibody Rapid Test Strip

Sample Type

Used to detect Toxoplasma antibodies in porcine serum, this test can be used as an adjunct in the diagnosis of porcine toxoplasmosis.

Product Introduction

The Porcine Toxoplasma Antibody Rapid Test Strip is based on the principle of the colloidal gold immunochromatographic assay (GICA). After the sample is added to the sample well, it migrates along the chromatographic membrane together with the colloidal gold conjugate. If anti-TOX antibodies are present in the sample, a purple-red color will appear on the test line, if no anti-TOX antibodies are present, no color reaction will occur.

Product Features

Delivers results quickly, making it ideal for rapid on-site screening.
Effectively improves the sensitivity and specificity of testing.
Requires no instruments or equipment, making it easier to use.

Instructions for Use

Avoid repeated freeze-thaw cycles of the samples.
Strictly follow the instructions in this manual.
Allow the kit to reach room temperature before use.

Kit Components

Test strips: 50
Instructions: 1 copy

Storage and Shelf Life

Store the kit in a dry environment at 2°C-30°C.
This product has a shelf life of 24 months. See the packaging for the production date.

Important Notes

Sample preparation is straightforward, and the procedure is simple, requiring minimal demands on the testing environment and personnel. Please follow the instructions or SOP for specific details.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY