

Impatiens Necrotic Spot Virus (INSV) RT-PCR Negative Control (10)



Cat. No. CPD-593

Lot. No. (See product label)

Product Name

Impatiens Necrotic Spot Virus (INSV) RT-PCR Negative Control (10)

Product Overview

Impatiens Necrotic Spot Virus (INSV) RT-PCR Negative Control (10) is designed for the detection of Impatiens Necrotic Spot Virus (INSV).

Kit Components

The kits contain a thermostable Reverse Transcriptase and a Hot-Start DNA Polymerase, pathogen specific primers and dNTPs, DNA-based positive and negative controls for the PCR reaction.

Detection method

RNA PCR

Assay Protocol

1. Isolated Virus-RNA from infected plant;
2. Reverse transcription + cDNA Amplication;
3. Amplified cDNA;
4. Detection in Agarose gel electrophoresis.

Sample Type

Impatiens

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY