MUT-1820 arg7-2 mt-

Cat. No. ALS-02025

Lot. No. (See product label)

Subategory

Mutants

Description

The arg7-2 mutation was induced using ethyl methanesulfonate. This allele can complement with the original arg7, arg7-5, arg7-8 (aka arg2), and arg7-13. It cannot complement with arg7-3, arg7-4, arg7-7 or arg7-12. It seems to be stable in culture.

Species

Chlamydomonas

Locus

ARG7

Chromosome

1

Phenotype

Requires arginine

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Address: SUITE 209, 17 Ramsey Road, Shirley, NY 11967 E-mail: info@ Tel: 1-631-356-5138 Fax: 1-631-910-2166 https://www.lifeasible.com/algae